

Greenfix**Manufacturer: Diapath** Via Savoldini,71- 24057 MARTINENGO- BG- PH+39.0363.986.411 info@diapath.com

CODE	PACKAGING
P0113	1 liter
P0116	5 liters
P0117	10 liters
P0118	20 liters

CODE	CONTAINER	PACKAGING	BOTTLE VOLUME	FILLED AT
Q28012	Container made of polypropylene with pressure cap	1x48 pz	30 ml	20 ml
Q28023		1x35 pz	60 ml	30 ml
Q28101		1x36 pz	150 ml	90 ml
Q28111		1x36 pz	250 ml	150 ml

Presentation

Ethandial based fixing solution, universal fixative formalin substitute.

Composition

ethandial	CAS No. 107-22-2	EC No 203-474-9
ethanol	CAS No 64-17-5	EC No 200-578-6

Method

Manual user available for further informations.

Soak the sample in a fixing solution.

Volumetric relationship between sample/fixative: 1:20**Fixation time:** Like formalin**Storage, handling and stability**

Store at room temperature. The product, if opportunely stored and integrally packed, has validity wrote on label MSDS available for further informations. Professional use.